

Step Stools



STS121805



T0800
(optional transport cart)

Designed with safety in mind, our step stools are a tip resistant, heavy-duty design. Single step stools can connect together to create a large, secure platform. The optional roll around transport cart allows for easy transportation and storage of up to 5 single step stools.

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at sales@macmedical.com.

Standard Models

Model #	Overall Dimensions
STS121805 (single step)	18"W x 12"D x 6.25"H
STS263015 (two steps)	30"W x 18"D x 15"H
T0800 (optional transport cart)	20"W x 19"D x 35"H

*Note: Additional sizes available for specific requirements.
Please contact us for further assistance.*

STANDARD FEATURES

- Heavy-duty, stainless steel construction
- Tip resistant
- Corrugated, black rubber top surface on each stool
- Rubber floor tips
- Single step stools (STS121805) are stackable and connect together to create a large platform

OPTIONAL FEATURES

- Easy roll around transport cart with 4 swiveling casters (2 locking) available for STS121805 – holds up to 5 single step stools



Corporate Office
325 West Main Street
Belleville, IL 62220

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

Phone: 618-476-3550
Toll Free: 877-828-9975
Fax: 618-476-3337
sales@macmedical.com
www.macmedical.com

MAC Medical, Inc.
Publication No. FFU-03 Rev B
October 2016
Information regarding this product is subject to change without prior notice.

MADE IN THE USA