

Prep & Pack Workstations

Ergonomically
Friendly



TWS-3072
(shown with multiple optional features)

STANDARD FEATURES

- 100% stainless steel construction
- Modular design allows for multiple configurations
- Integrated footrest
- 4" casters (2 swivel & 2 locking)

OPTIONAL FEATURES

- **Ergonomically friendly** adjustable heights (manual or electronic)
 - Manual adjustable height range: 32"–43" (working surface)
 - Electronic adjustable height range: 32"–42" (working surface)
- Bin rail
- Overhead light
- Articulating magnification light
- Electrical strip
- CPU holder
- Pull out keyboard tray
- Magnetic dry erase board
- Adjustable flat shelf for riser frame
- Adjustable angled shelf for riser frame
- Lower shelf
- Side mount waste basket holder with waste basket
- Spool holder
- Drawer boxes
 - Single 6" drawer box
 - Double 6" drawer box
 - Double drawer box with 6" drawer and file drawer

Designed with ergonomics in mind, the Prep & Pack Workstations feature a modular design, allowing for multiple configurations and adjustability for all body types. Choose from the many optional features including adjustable height (manual or electronic), overhead light, magnification light, CPU holder, electrical strip, and more.

Standard Models

Base Model #	Overall Dimensions (height at lowest setting)
TWS-3060	60"W x 32.25"D x 81"H
TWS-3072	72"W x 32.25"D x 81"H
TWS-3660	60"W x 38.25"D x 81"H
TWS-3672	72"W x 38.25"D x 81"H

Note: Additional sizes available for specific requirements. Please contact us for further assistance.

Note: Model # will change based on features added.



Corporate Office
325 West Main Street
Belleville, IL 62220

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

Phone: 618-476-3550
Toll Free: 877-828-9975
Fax: 618-476-3337
sales@macmedical.com
www.macmedical.com

MAC Medical, Inc.
Publication No. FFU-06 Rev C
October 2016
Information regarding this product is subject to change without prior notice.

MADE IN THE USA