

*New*

# Cool Touch Handles

Eliminate "hot handle" discomfort from cart washing



**MMCV95**  
Vertical Handle Case Cart  
(with optional Cool Touch Handles)



Eliminate the discomfort of "hot handles" from cart washers by using our new safe and comfortable Cool Touch Handles. These optional handles are designed for use on any of our Vertical Handle Case Carts (or any 1" diameter vertical handle). The retrofit kit includes two Cool Touch Handles with corresponding hardware.


## STANDARD FEATURES

- Comfortable & easy grip
- Made of durable Santoprene™ TPV material
- Simple assembly (hardware included with each kit)
- High-temperature sustainability
- 5.625"H x 2" diameter

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at **sales@macmedical.com**.

### COOL TOUCH HANDLE KIT (Part # WC011-100-KT)

Kit includes two Cool Touch Handles with hardware



**Corporate Office**  
325 West Main Street  
Belleville, IL 62220

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550  
Toll Free: 877-828-9975  
Fax: 618-476-3337  
sales@macmedical.com  
www.macmedical.com

MAC Medical, Inc.  
Publication No. FFU-17 Rev A  
December 2016  
Information regarding this product is  
subject to change without prior notice.