

Scope Cabinets

Efficient drying and storage for multiple endoscopes



IC184292
(shown with optional power vent)

STANDARD FEATURES

- Fully welded construction
- Glass hinged doors with louvered vents
- Scope hangers
- Lower rack to keep scopes aligned
- Removable drain pan
- Padded back wall to prevent damage

OPTIONAL FEATURES


- Power vent (includes ventilation fan, digital time switch, and HEPA filter)
- Sloped top (adds 6" to total height)

Standard Models

Base Model #	Overall Dimensions	# of Scope Hangers
IC183692	18"D x 36"W x 92"H	9 or 18
IC184292	18"D x 42"W x 92"H	11 or 22

Note: Additional sizes available for specific requirements, please contact us for further assistance

Note: Model # will change based on features added



MADE USA
IN THE

Customer needs are our first priority.

mac
MEDICAL

Corporate Office
325 West Main Street
Belleville, IL 62220

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at sales@macmedical.com.

Phone: 618-476-3550
Toll Free: 877-828-9975
Fax: 618-476-3337
sales@macmedical.com
www.macmedical.com

MAC Medical, Inc.
Publication No. FFU-18 Rev A
August 2017
Information regarding this product is subject to change without prior notice.