

# General Transport Stretcher PT-1000 is now PT1001

## Cross Reference Guide



**550**  
POUND  
Weight  
Capacity

### PT-1000 *Old Model*

#### STANDARD FEATURES

- 550 lb. weight capacity
- 4 brake/steer pedals
- Single Fowler assisted back rest with 90° adjustability
- 4" standard mattress
- Storage area underneath for patient belongings
- 8" Tente casters (1 directional & 3 swiveling)
- Height adjustment range: 22"-35" (bottom of mattress)
- 26", 29", and 31" bed width (PT-1000, PT-1000-29, and PT-1000-31)
- Dual pedestal lift mechanism
- Lift/Trendelenburg pedals on both sides
- Trendelenburg & reverse Trendelenburg +/- 18°
- Fold forward side rails with built in push handles
- Full perimeter bumper
- 2 hook removable telescopic IV pole
- 6 IV pole sockets available

#### OPTIONAL FEATURES

- Envelope style rail bumper pads
- Fold down push handles
- Detachable foot board
- Paper roll dispenser kit
- Knee gatch

*Redesigned*



**750**  
POUND  
Weight  
Capacity

### PT1001 *New Model*

#### STANDARD FEATURES

- 750 lb. weight capacity *New!*
- 6 brake/steer pedals *New!*
- Dual Fowler assisted lift back rest with 90° adjustability *New!*
- 4" comfort mattress *New!*
- Larger storage area underneath for patient belongings *New!*
- Ergonomically located oxygen tank holder *New!*
- 4 restraint belt locations *New!*
- 8" Tente casters (1 directional & 3 swiveling)
- Height adjustment range: 22"-35" (bottom of mattress)
- 29" or 31" bed width (PT1001 & PT1001-31)
- Dual pedestal lift mechanism
- Lift/Trendelenburg pedals on both sides
- Trendelenburg & reverse Trendelenburg +/- 18°
- Fold forward side rails with built in push handles
- Full perimeter bumper
- 2 hook removable telescopic IV pole
- 6 IV pole sockets available

#### OPTIONAL FEATURES

- 5th wheel add-on kit *New!*
- Envelope style rail bumper pads
- Fold down push handles
- Detachable foot board
- Paper roll dispenser kit
- Knee gatch

*Customer needs are our first priority.*

MAC Medical, Inc.  
Printed in USA  
Publication No. CRFU-01 Rev B  
October 2016

Information regarding this product is  
subject to change without prior notice.



**Corporate Office**  
325 West Main Street  
Belleville, IL 62220

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • Fax: 618-476-3337  
Email: sales@macmedical.com • Website: www.macmedical.com