

*The next generation of warming cabinets*

# D-Series Data Logging Blanket & Fluid Warming Cabinets



The D-Series (Data Logging) Blanket & Fluid Warming Cabinets provide independent, digitally controlled heating chambers that offer actual temperature and set point displays. Data is recorded every 30 minutes and stored for up to two years. Choose from single, dual, or triple chamber units.

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at **sales@macmedical.com**.

## STANDARD FEATURES

- "Easy Log" integrated data logger with removable USB – *now standard!*
- Keyed door lock for single & upper chamber units – *now standard!*
- Adjustable feet – *now standard!*
- Automatically records data every 30 minutes
- Internal storage up to 2 years
- Displays actual temperature and set points at all times
- Independent digitally controlled heating chambers
- Rapid warm time (2-6 hours depending on load)
- Keyed temperature lock out
- Push button digital operation
- Large digital display(s)
- Fully insulated to provide uniform heating
- Audible & visual over-temp alarm

## OPTIONAL FEATURES

- Full view glass door
- Left hand door hinge
- Mobile base or stand
- 220/240V power supply
- Celsius
- Two (2) recessed unit options
- Seismic braces
- Roll out basket(s) with cool touch handles
- Pass through chamber
- Sloped top
- IV/Injection fluids
- Keyed door lock for intermediate & lower chamber units
- Electronic keypad lock
- Custom sizes available



**MADE IN THE USA**

[www.macmedical.com](http://www.macmedical.com)


# How to Order D-Series

Model #	Overall Size	Chambers	Door Type	Door Hinge	Base Style	Other Options
SWC151822	15.5"D x 18"W x 22"H	Single	G = Glass Door	LH = Left Hinge	† 2B = 2" Base	220 = 220/240V Power Supply
SWC182424	20.5"D x 24"W x 24.5"H	Single	Blank = Stainless Steel Door (Standard)	Blank = Right Hinge (Standard)	† 4B = 4" Base	C = Celsius
SWC242424	26.5"D x 24"W x 24.5"H	Single			† 6B = 6" Base	R1 = Recessed Unit with Insulation Wrap, No Top Panel, No Side Panels, and No Trim Kit
SWC183024	20.5"D x 30"W x 24.5"H	Single			NB = No Base	R2 = Recessed Unit with Top Panel, Side Panels, and Trim Kit
SWC243024	26.5"D x 30"W x 24.5"H	Single			MB = Mobile Base	SB = Seismic Braces
SWC182436	20.5"D x 24"W x 36" H	Single			WB = Roll Out Basket	P = Pass Through Chamber
SWC183036	20.5"D x 30"W x 36"H	Single			IV = IV/Injection Fluids	† ST = Sloped Top
SWC242436	26.5"D x 24"W x 36"H	Single			†† IDL = Intermediate Chamber Door Lock	††† LDL = Lower Chamber Door Lock
SWC243036	26.5"D x 30"W x 36"H	Single			†††† EL = Electronic Keypad Lock	
SWC182464	20.5"D x 24"W x 64.75"H	Single				
SWC242464	26.5"D x 24"W x 64.75"H	Single				
SWC183064	20.5"D x 30"W x 64.75"H	Single				
SWC243064	26.5"D x 30"W x 64.75"H	Single				
SWC182474	20.5"D x 24"W x 74.5"H	Single				
SWC183074	20.5"D x 30"W x 74.5"H	Single				
SWC242474	26.5"D x 24"W x 74.5"H	Single				
SWC243074	26.5"D x 30"W x 74.5"H	Single				
DWC182464T*	20.5"D x 24"W x 64.75"H	Dual				
DWC242464T*	26.5"D x 24"W x 64.75"H	Dual				
DWC183064T*	20.5"D x 30"W x 64.75"H	Dual				
DWC243064T*	26.5"D x 30"W x 64.75"H	Dual				
DWC182474T*	20.5"D x 24"W x 74.5"H	Dual				
DWC242474T*	26.5"D x 24"W x 74.5"H	Dual				
DWC183074T*	20.5"D x 30"W x 74.5"H	Dual				
DWC243074T*	26.5"D x 30"W x 74.5"H	Dual				
DWC182474E**	20.5"D x 24"W x 74.5"H	Dual				
DWC183074E**	20.5"D x 30"W x 74.5"H	Dual				
DWC242474E**	26.5"D x 24"W x 74.5"H	Dual				
DWC243074E**	26.5"D x 30"W x 74.5"H	Dual				
TWC183078	20.5"D x 30"W x 78.75"H	Triple				
TWC243078	26.5"D x 30"W x 78.75"H	Triple				

the **NEXT** generation of warming cabinets

† Not available for SWC151822  
 †† Only available for Triple Chamber Units  
 ††† Only available for Dual and Triple Chamber Units  
 †††† Only available for Single Chamber Units

\* T = Chambers in Thirds      \*\* E = Equal Chambers



**Example 1:** Model # SWC182464, Glass Door, Left Hinge, Mobile Base, Sloped Top = Order Model # SWC182464-G-LH-MB-ST

**Example 2:** Model # DWC243074E, Stainless Steel Door, Right Hinge, 4" Base, Celsius, Lower Chamber Door Lock = Order Model # DWC243074E-4B-C-LDL

Note: Additional sizes available for specific requirements. Please contact us for further assistance.



Customer needs are our first priority.

MAC Medical, Inc.  
 Printed in USA  
 Publication No. FFU-02 Rev F  
 January 2018

Information regarding products are subject to change without prior notice.



Corporate Office  
 325 West Main Street  
 Belleville, IL 62220

Manufacturing Plant  
 820 South Mulberry Street  
 Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • www.macmedical.com