
Surgical Scrub Sinks

ES-Series



Model # ES47 (Dual Station)

DIMENSIONS



NOTES:

- 1. MATERIAL: 14 GAUGE 304 STAINLESS STEEL
- 2. ALL PLUMBING 1/2 COPPER OR BRASS PIPE.
- 3. UTILITY REQUIREMENT

 COLD WATER 20 TO 50 PSIG 70 DEG F MAX

 HOT WATER 20 TO 50 PSIG 120 DEG F TO 140 DEG F
- 4. POWER REQUIREMENTS (FOR SENSOR OPERATIONS) 120 VAC, 60 Hz, SINGLE PHASE, 3.0 AMP OR 220 VAC, 60 Hz, SINGLE PHASE, 1.5 AMP. INSTALLED PER LOCAL BUILDING CODES.
- 5. NET WEIGHT OF SINK 110 LB.

Surgical Scrub Sinks

ES-Series



Model # ES47 (Dual Station)
with Wall Mounting Brackets (included)

DIMENSIONS

