

# Processing Sinks

Pipe Base Units

**Ergonomically  
Friendly**



**PS2612035-9**  
(electronic adjustable height shown with left to right workflow and optional pre-rinse faucet, DI water, air gun, and water gun)



Note: ETL compliance for adjustable height models only

## STANDARD FEATURES

- Open pipe base
- Single, dual, or triple basins
- Fully welded construction with reinforced tops
- Right to left workflow or left to right workflow
- 37" standard height for stationary sinks
- 4" backsplash
- Under shelf
- Adjustable feet
- Flat working surface
- Etched gallon markings
- Basin size: 24"W x 17"L x 10"D
- Full CAD rendering of specified sink

## OPTIONAL FEATURES

- **Ergonomically friendly** electronic adjustable height with four position programming
  - Electronic adjustable height range: 35"-42"
- Wrap around apron
- Mixing valve faucet
- Pre-rinse faucet
- Deionized (DI) water faucet
- Air gun kit (with 5 tips)
- Water gun kit (with 8 tips)
- Shut off valve for air/water guns
- Eye wash
- Drain boards
- Larger basin: 30"W x 17"L x 10"D
- Backsplash (6", 8" & 12")
- Side splash (4", 6", 8" & 12")
- Tubing manifold
- Hand lever drains
- Foot operated drains (not available for adjustable height units)
- Copper drain system (not compatible with deionized water)
- Overhead LED hood

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at **sales@macmedical.com**.

## Standard Pipe Base Models (Stationary)

| Base Model # | Overall Dimensions                |
|--------------|-----------------------------------|
| PS267237-1   | 26"D x 72"W x 37"H, Single Basin  |
| PS269637-1   | 26"D x 96"W x 37"H, Dual Basin    |
| PS2612037-1  | 26"D x 120"W x 37"H, Triple Basin |

Note: Additional sizes available for specific requirements, please contact us for further assistance

Note: Model # will change based on features added

## Standard Pipe Base Models (Electronic Adjustable Height)

| Base Model # | Overall Dimensions                |
|--------------|-----------------------------------|
| PS267235-9   | 26"D x 72"W x 35"H, Single Basin  |
| PS269635-9   | 26"D x 96"W x 35"H, Dual Basin    |
| PS2612035-9  | 26"D x 120"W x 35"H, Triple Basin |

Note: Additional sizes available for specific requirements, please contact us for further assistance

Note: Model # will change based on features added



**Corporate Office**  
325 West Main Street  
Belleville, IL 62220

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550  
Toll Free: 877-828-9975  
Fax: 618-476-3337  
sales@macmedical.com  
www.macmedical.com

MAC Medical, Inc.  
Publication No. FFU-20 Rev A  
April 2018  
Information regarding this product is subject to change without prior notice.

**MADE IN THE USA**