


# 5th Wheel Transport Stretcher PTF1000



**750**  
POUND  
Weight  
Capacity

Retractable  
Fifth Wheel

The PTF1000 5th Wheel Transport Stretcher allows for easy and comfortable patient transfer with a weight capacity of up to 750 lbs. The retractable center fifth wheel combines effortless directional movement and cornering with 4 swiveling wheels for ease of use, control, and safety.

## Contact Us

Contact our customer service department toll free at 877-828-9975, or by email at [sales@macmedical.com](mailto:sales@macmedical.com)

**MADE  
IN THE USA**

## STANDARD FEATURES

- Retractable center fifth wheel provides effortless directional movement and cornering
- 750 lb. weight capacity
- Height adjustment range: 22"-35" (bottom of mattress)
- Ergonomically located oxygen tank holder
- Effortless steering, pivoting, and cornering with 8" Tente casters (4 swiveling)
- 6 brake/steer pedals
- Dual pedestal lift mechanism
- Lift/Trendelenburg pedals on both sides
- Trendelenburg & reverse Trendelenburg +/- 18°
- Dual Fowler assisted lift back rest with 90° adjustability
- Fold forward side rails with built in push handles
- 29" or 31" bed width
- 4" comfort mattress
- 2 hook removable telescopic IV pole
- 6 IV pole sockets available
- Large storage area underneath for patient belongings
- Full perimeter bumper to protect your walls and fixtures
- 4 restraint belt locations

## OPTIONAL FEATURES

- Fold down push handles
- Knee gatch
- Paper roll dispenser kit \*
- Detachable foot board \*
- Envelope style rail bumper pads \*

\* Additional lead time required (7-10 business days)

Model # will change based on features added

## PRODUCT SPECIFICATIONS

### PTF1000

- 29" bed width
- Overall dimensions: 83 5/16" long x 33 9/16" wide
- Patient surface dimensions: 75 7/16" long x 29 1/8" wide

### PTF1001-31

- 31" bed width
- Overall dimensions: 83 5/16" long x 35 9/16" wide
- Patient surface dimensions: 75 7/16" long x 31 1/8" wide



[www.macmedical.com](http://www.macmedical.com)

## Features




Retractable center fifth wheel provides effortless directional movement and cornering

Dual Fowler assisted lift back rest with 90° adjustability



Ergonomically located oxygen tank holder




Large storage area for patient belongings




6 brake/steer pedals (two on each side & one on each head and foot end)



Effortless steering, pivoting, and cornering with 8" Tente casters (4 swiveling)




Trendelenburg & reverse Trendelenburg +/- 18°




Height adjustment range: 22"-35" (bottom of mattress)

## Dimensions


### PTF1000



### PTF1000-31



### Height Adjustment Ranges (Both Models)



INSPIRED  
TO BE BETTER

Customer needs are our first priority.

MAC Medical, Inc.  
Printed in USA  
Publication No. FFU-01 Rev D  
July 2021

Information regarding this product is  
subject to change without prior notice.



**Corporate Office**  
200 Admiral Weinel Blvd.  
Columbia, Illinois 62236

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • [www.macmedical.com](http://www.macmedical.com)