

# Step Stools



STS121805



T0800  
(optional transport cart)

Designed with safety in mind, our step stools are a tip resistant, heavy-duty design. Single step stools can connect together to create a large, secure platform. The optional roll around transport cart allows for easy transportation and storage of up to 5 single step stools.

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at [sales@macmedical.com](mailto:sales@macmedical.com).

## Standard Models

Model #	Overall Dimensions
STS121805 (single step)	18"W x 12"D x 6.25"H
STS263015 (two steps)	30"W x 18"D x 15"H
T0800 (optional transport cart)	20"W x 19"D x 35"H

*Note: Additional sizes available for specific requirements.  
Please contact us for further assistance.*

## STANDARD FEATURES

- Heavy-duty, stainless steel construction
- Tip resistant
- Corrugated, black rubber top surface on each stool
- Rubber floor tips
- Single step stools (STS121805) are stackable and connect together to create a large platform

## OPTIONAL FEATURES

- Easy roll around transport cart with 4 swiveling casters (2 locking) available for STS121805 – holds up to 5 single step stools



**MADE**  
IN THE **USA**

**Corporate Office**  
200 Admiral Weinel Blvd.  
Columbia, Illinois 62236

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550  
Toll Free: 877-828-9975  
Fax: 618-476-3337  
[sales@macmedical.com](mailto:sales@macmedical.com)  
[www.macmedical.com](http://www.macmedical.com)

MAC Medical, Inc.  
Publication No. FFU-03 Rev C  
July 2021  
Information regarding this product is subject to change without prior notice.