Premium upgrade with smart technology TS-Series Touch Screen **Blanket & Fluid Warming Cabinets**





The TS-Series (Touch Screen) Blanket & Fluid Warming Cabinets provide the perfect blend of smart technology, value, and customer support in the industry. The intuitive touch screen gives you complete system control of your warmer, while the Ethernet/Wi-Fi capability allows remote connectivity (via smart phone app) giving you access to your warmer wherever you are.

Contact Us

Contact our customer service department toll free at 877-828-9975, or by email at sales@macmedical.com







STANDARD FEATURES

Many of the same features as the D-Series, but also includes:

- Intuitive touch screen control
- Intuitive software (supports multiple languages)
- FDA compliant data logging with digital signatures
- View live trends
- Alarm logging
- Multi-level user security
- Operator audit trail
- Built in Ethernet
- Wi-Fi capable (router to be supplied by end user)
- Remote monitoring & control (VNC)
- Remote tech support & diagnostics
- Text and email alarm notifications
- File transfer via USB or FTP web page
- Smart phone/device access

OPTIONAL FEATURES

- Full view glass door
- Left hand door hinge
- · Mobile base or stand
- 220/240V power supply
- Celsius
- Two (2) recessed unit options
- Seismic braces
- Roll out basket(s) with cool touch handles
- Pass through chamber
- Sloping top
- Keyed door lock for lower chamber units
- Electronic keypad lock















Check on the status of your warmer(s) with your smart phone device



Manage multiple warmers remotely or locally, via Ethernet or Wi-Fi (router for Wi-Fi to be supplied by end user)




Optional roll out basket makes organization and item selection easy. Cool touch handles keep your hands from getting burned.



Base Model #	Overall Size	Chambers	Door Type	Door Hinge	Base Style	Other Options
SWC182424-TS	20.5"D x 24"W x 24.5"H	Single	G = Glass Door Blank = Stainless Steel Door (Standard) LH = Left Hinge Blank = Right Hinge (Standard)	Blank = Right	2B = 2" Base	220 = 220/240V Power Supply
SWC243024-TS	26.5"D x 30"W x 24.5"H	Single			4B = 4" Base	
SWC243036-TS	26.5"D x 30"W x 36"H	Single				C = Celsius
SWC182474-TS	20.5"D x 24"W x 74.5"H	Single			6B = 6" Base	
SWC243074-TS	26.5"D x 30"W x 74.5"H	Single		05 0 5030	R1 = Recessed Unit	
DWC183074T-TS*	20.5"D x 30"W x 74.5"H	Dual			NB = No Base	with Insulation Wrap No Top Panel, No Sic Panels, and No Trim
DWC243074T-TS*	26.5"D x 30"W x 74.5"H	Dual				
DWC243074E-TS**	26.5"D x 30"W x 74.5"H	Dual			MB = Mobile Base	R2 = Recessed Unit with Top Panel, Side Panels, and Trim Kit SB = Seismic Braces WB = Roll Out Basket P = Pass Through Chamber
			S	MAR1 technol	logy	EL = Electronic Keypad Lock ST = Sloping Top
						LDL = Lower Chamber Door Lock (Only available for Dual Chamber Units)

* T = Chambers in Thirds



Example 1: Model # DWC243074E-TS, Stainless Steel Door, Right Hinge, 4" Base, Celsius, Lower Chamber Door Lock = Order Model # DWC243074E-TS-4B-C-LDL

Example 2: Model # SWC182474-TS, Glass Door, Left Hinge, Sloping Top, Mobile Base = Order Model # SWC182474-TS-G-LH-ST-MB

Note: Model # will change based on features added











Customer needs are our first priority.



Corporate Office 200 Admiral Weinel Blvd. Columbia, Illinois 62236

Manufacturing Plant 820 South Mulberry Street Millstadt, IL 62260