

Prep & Pack Workstations

Ergonomical and versatile workstation for all body types

**Ergonomically
Friendly**

Designed with ergonomics in mind, the Prep & Pack Workstations feature a modular design, allowing for multiple configurations and adjustability for all body types. Choose from the many optional features including adjustable height (manual or electronic), overhead light, articulating magnification light, CPU holder, electrical strip, and more.

STANDARD FEATURES

- 100% stainless steel construction
- Modular design allows for multiple configurations
- Integrated footrest
- 5" casters (2 swivel & 2 locking)

OPTIONAL FEATURES

- **Ergonomically friendly** adjustable heights (manual or electronic)
 - Adjustable height range: 36"–48"
- Bin rail
- Overhead light
- Articulating magnification light
- Electrical strip
- CPU holder
- Pull out keyboard tray
- Magnetic dry erase board
- Adjustable flat shelf for riser frame
- Adjustable angled shelf for riser frame
- Lower shelf
- Spool holder
- Drawer boxes (single or double)



TWS-3672-EAH-B
(shown with multiple optional features)

Standard Models

Base Model #	Overall Dimensions (height at lowest setting)
TWS-3060-B	60"W x 32"D x 84"H
TWS-3072-B	72"W x 32"D x 84"H
TWS-3660-B	60"W x 38"D x 84"H
TWS-3672-B	72"W x 38"D x 84"H

Note: Model # will change based on features added.

MADE IN THE USA

mac
MEDICAL

Corporate Office
200 Admiral Weinle Blvd.
Columbia, Illinois 62236

Phone: 618-476-3550
Toll Free: 877-828-9975
Fax: 618-476-3337
sales@macmedical.com
www.macmedical.com

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

MAC Medical, Inc.
Publication No. FFU-06 Rev I
March 2024
Information regarding this product is
subject to change without prior notice.