

# Stainless Steel Case Carts



MAC Medical, Inc. manufactures a full line of stainless steel case carts that feature fully welded construction for strength and durability on all carts, as well as easy service bolt-on features on closed case carts. Choose from several designs including Closed, Sealed, and Open Case Carts. Our engineering capabilities allow us to custom design any specific case cart to satisfy any need your facility may have.

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at **[sales@macmedical.com](mailto:sales@macmedical.com)**.

[Closed Case Carts](#)

---

[Sealed Case Carts](#)

---

[Open Case Carts](#)

---

[View more information on the back >>](#)

**MADE**  
IN THE **USA**



[www.macmedical.com](http://www.macmedical.com)



## Closed Case Carts

### STANDARD FEATURES

- Fully welded construction with inner skeletal system for strength and our easy service bolt-on features
- Soundproof tops and doors for quiet operation
- 6" maintenance-free, high temp, high pressure stainless steel casters with sealed ball bearings (3 swivel & 1 fixed directional)
- Bolt-on exterior push handle for better control
- Bolt-on full wrap-around non-marring bumper assembly
- Bolt-on door hinges
- Bolt-on positive lock door latches
- Bolt-on 8.5" x 11" card holder
- Bolt-on tamper evidence security hasps

### OPTIONAL FEATURES

- Glass door(s)
- Etching for serial tag ID/tracking
- Etching for end user logo
- Adjustable solid shelf (non roll out)
- Adjustable solid shelf (roll out)
- Adjustable perforated shelf (non roll out)
- Adjustable perforated shelf (roll out)
- Adjustable wire shelf (non roll out)
- Adjustable wire shelf (roll out)
- Side extensions
- Flag system (one each of red, yellow and green)
- Clean/dirty indicator (slide or wheel)
- Total locking casters
- Tow hitch



## Sealed Case Carts

### STANDARD FEATURES

- Fully welded construction
- Soundproof tops and doors for quiet operation
- Exterior push handle for better control
- Door gaskets for increased dust control
- Locking rods to keep doors secure
- Full wrap around non-marring heavy-duty bumper
- Magnetic door catches to keep doors open while in use
- 8.5" x 11" card holder
- Tamper evidence security tabs
- 6" maintenance-free, high-temperature, high-pressure stainless steel casters with sealed ball bearings (2 rigid & 2 swivel with brake)

### OPTIONAL FEATURES

- Etching for serial tag ID/tracking
- Etching for end user logo
- Adjustable solid shelf (non roll out)
- Adjustable solid shelf (roll out)
- Adjustable perforated shelf (non roll out)
- Adjustable perforated shelf (roll out)
- Adjustable wire shelf (non roll out)
- Adjustable wire shelf (roll out)
- Side extensions
- Flag system (one each of red, yellow and green)
- Clean/dirty indicator (slide or wheel)
- Total locking casters
- Tow hitch



## Open Case Carts

### STANDARD FEATURES

- Open design
- Sturdy built cart – perfect for storage, supplies, linens, or distribution
- Full wrap around non-marring bumper
- Stainless steel guard rail
- 6" composite directional locking casters

### OPTIONAL FEATURES

- Adjustable solid shelf (non roll out)
- Adjustable perforated shelf (non roll out)
- Adjustable wire shelf (non roll out)
- Fixed shelving
- Wire grid on back side
- 6" maintenance-free, high temp, high pressure stainless steel casters with sealed ball bearings (3 swivel & 1 fixed directional)
- Total locking casters
- 6" plastic bin
- 10" plastic bin

For questions, or to place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at [sales@macmedical.com](mailto:sales@macmedical.com).

INSPIRED  
TO BE BETTER

*Customer needs are our first priority.*

MAC Medical, Inc.  
Printed in USA  
Publication No. FFU-07 Rev G  
July 2022

Information regarding this product is subject to change without prior notice.



**Corporate Office**  
200 Admiral Weinel Blvd.  
Columbia, Illinois 62236

**Manufacturing Plant**  
820 South Mulberry Street  
Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • [www.macmedical.com](http://www.macmedical.com)