


Surgical Scrub Sinks

SS-Series



OSHDP
Pre-Approved

The SS-Series is our premier line of surgical scrub sinks. Available in single, dual, and triple basins, these sinks combine the best material, components, and craftsmanship available on the market today.

STANDARD FEATURES

- Single, dual, or triple stations
- Deep sloping basin
- Adjustable thermostatic controlled pressure regulating mixing valve (one per station)
- Knee operated water & soap (optional with IR sinks)
- Check stop strainers for ease of maintenance
- High mount swivel gooseneck spout with adjustable swivel spray aerator
- Flow control valve
- Backflow preventers
- Infrared sinks are equipped with a 24 volt plug in transformer
- Splash screen dividers with double and triple station sinks
- Main support and Z-brackets (wall support optional)
- 304 stainless steel polished and ground to a #4 finish

OPTIONAL FEATURES

- Infrared operation
- Infrared soap dispenser
- Manual bypass to add knee operated water & soap with IR operation
- Knee operated soap with IR water operation
- Foot operated water
- Foot operated soap
- Digital scrub timer
- Eye wash
- Thermostatic mixing valve (add on to eye wash option)
- Laminar flow
- 220/240V power supply
- Pedestal base
- In-wall chair carrier
- Wall shelf over sink
- Lighted wall shelf over sink



SS32
Single Station



SS64
Dual Station



SS96
Triple Station

For more information, or to place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at sales@macmedical.com.



MADE
IN THE **USA**

www.macmedical.com

SS-Series Features



OSHDP
Pre-Approved

SS96
(shown with multiple optional features)

How to Order SS-Series

Infrared Operation	Infrared Soap Dispenser	Manual Bypass to add Knee Operated Water and Soap with IR Operation	Knee Operated Water with IR Operation	Knee Operated Soap with IR Operation	Foot Operated Water	Foot Operated Soap	Digital Scrub Timer	Eye Wash	Thermostatic Mixing Valve (add on to Eye Wash option)	Laminar Flow	220/240V	Pedestal Base
IR	IRD	MB	KOW	KSP	FOW	FSP	TM	EYE	TV	LF	220	PB

Note: Model # will change based on features added

Example 1: Base Model # SS32 with Infrared Soap Dispenser, Digital Scrub Timer, and Pedestal Base = Model # **SS32-IRD-TM-PB**

Example 2: Base Model # SS96 with Foot Operated Water, Eye Wash, and Laminar Flow = Model # **SS96-FOW-EYE-LF**

Recommended Optional Products

Model #	In-Wall Chair Carrier #	Shelf Over Sink #	Lighted Shelf Over Sink #
SS32	S0001	SH0001 (11.5" D x 32" W)	SH0001-L
SS64	S0002	SH0002 (11.5" D x 64" W)	SH0002-L
SS96	S0003	SH0003 (11.5" D x 32" W)	SH0003-L



Optional Infrared Soap Dispenser

Standard SS-Series Models

Base Model #	Overall Dimensions
SS32	27.5"D x 32"W x 37"H
SS64	27.5"D x 64"W x 37"H
SS96	27.5"D x 96"H x 37"H

Note: Additional sizes available for specific requirements, please contact us for further assistance

Note: Model # will change based on features added

INSPIRED
TO BE BETTER

Customer needs are our first priority.

MAC Medical, Inc.
Printed in USA
Publication No. FFU-14 Rev C
August 2021



Corporate Office
200 Admiral Weinel Blvd.
Columbia, Illinois 62236

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • Fax: 618-476-3337
Email: sales@macmedical.com • Website: www.macmedical.com

Information regarding this product is subject to change without prior notice.