

General Transport Stretcher PT1001

750
POUND
Weight
Capacity



STANDARD FEATURES

- 750 lb. weight capacity
- Height adjustment range: 22"-35" (bottom of mattress)
- Ergonomically located oxygen tank holder
- Effortless steering, pivoting, and cornering with 8" Tente casters (1 directional & 3 swiveling)
- 6 brake/steer pedals
- Dual pedestal lift mechanism
- Lift/Trendelenburg pedals on both sides
- Trendelenburg & reverse Trendelenburg +/- 18°
- Dual Fowler assisted lift back rest with 90° adjustability
- Fold forward side rails with built in push handles
- 29" or 31" bed width
- 4" comfort mattress
- 2 hook removable telescopic IV pole
- 6 IV pole sockets available
- Large storage area underneath for patient belongings
- Full perimeter bumper to protect your walls and fixtures
- 4 restraint belt locations

OPTIONAL FEATURES

- 5th wheel add-on kit
- Fold down push handles
- Knee gatch
- Paper roll dispenser kit *
- Detachable foot board *
- Envelope style rail bumper pads *

* Additional lead time required (7-10 business days)

PRODUCT SPECIFICATIONS

PT1001

- 29" bed width
- Overall dimensions: 83 5/16" long x 33 9/16" wide
- Patient surface dimensions: 75 7/16" long x 29 1/8" wide

PT1001-31

- 31" bed width
- Overall dimensions: 83 5/16" long x 35 9/16" wide
- Patient surface dimensions: 75 7/16" long x 31 1/8" wide

The PT1001 allows for easy and comfortable patient transfer with a weight capacity of up to 750 lbs. Economical and versatile, the PT1001 is an ideal solution for any department where cost efficient quality is required.


For questions or to place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at sales@macmedical.com.

MADE IN THE USA

mac
MEDICAL

www.macmedical.com

Features



Effortless steering, pivoting, and cornering with 8" Tente casters (1 directional & 3 swiveling)


Dual Fowler assisted lift back rest with 90° adjustability




Large storage area for patient belongings




Ergonomically located oxygen tank holder




6 brake/steer pedals (two on each side & one on each head and foot end)



Trendelenburg & reverse Trendelenburg +/- 18°




Height adjustment range: 22"-35" (bottom of mattress)




5th wheel add-on kit (optional)

Dimensions


PT1001



PT1001-31



Height Adjustment Ranges (Both Models)



INSPIRED
TO BE BETTER

Customer needs are our first priority.

MAC Medical, Inc.
 Printed in USA
 Publication No. FFU-16 Rev C
 August 2021



Corporate Office
 200 Admiral Weibel Blvd.
 Columbia, Illinois 62236

Manufacturing Plant
 820 South Mulberry Street
 Millstadt, IL 62260

Information regarding this product is subject to change without prior notice.

Phone: 618-476-3550 • Toll Free: 877-828-9975 • Fax: 618-476-3337
 Email: sales@macmedical.com • Website: www.macmedical.com