

Cool Touch Handles

Eliminate "hot handle" discomfort from cart washing



MMCV95
Vertical Handle Case Cart
(with optional Cool Touch Handles)



Eliminate the discomfort of "hot handles" from cart washers by using our new safe and comfortable Cool Touch Handles. These optional handles are designed for use on any of our Vertical Handle Case Carts (or any 1" diameter vertical handle). The retrofit kit includes two Cool Touch Handles with corresponding hardware.


STANDARD FEATURES

- Comfortable & easy grip
- Made of durable Santoprene™ TPV material
- Simple assembly (hardware included with each kit)
- High-temperature sustainability
- 5.625"H x 2" diameter

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at sales@macmedical.com.

COOL TOUCH HANDLE KIT (Part # WC011-100-KT)

Kit includes two Cool Touch Handles with hardware



Corporate Office
200 Admiral Weinel Blvd.
Columbia, Illinois 62236

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

Phone: 618-476-3550
Toll Free: 877-828-9975
Fax: 618-476-3337
sales@macmedical.com
www.macmedical.com

MAC Medical, Inc.
Publication No. FFU-17 Rev B
August 2021
Information regarding this product is
subject to change without prior notice.