

Scope Cabinets

Efficient drying and storage for multiple endoscopes



MSC-A2-4218
(shown with optional power vent)

STANDARD FEATURES

- Fully welded construction
- Glass hinged doors with louvered vents
- Scope hangers
- Lower rack to keep scopes aligned
- Removable drain pan
- Padded back wall to prevent damage

OPTIONAL FEATURES


- Power vent (includes ventilation fan, digital time switch, and HEPA filter)
- Sloping top (adds 6" to total height)

Standard Models

Base Model #	Overall Dimensions	# of Scope Hangers
MSC-A2-3618	18"D x 36"W x 92"H	9
MSC-A2-4218	18"D x 42"W x 92"H	11
MSC-A2-3624	24"D x 36"W x 92"H	18
MSC-A2-4224	24"D x 42"W x 92"H	22

Additional sizes available for specific requirements, please contact us for further assistance

Model # will change based on features added



Contact Us

Contact our customer service department toll free at **877-828-9975**, or by email at sales@macmedical.com



Corporate Office
200 Admiral Weinel Blvd.
Columbia, Illinois 62236

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

Phone: 618-476-3550
Toll Free: 877-828-9975
Fax: 618-476-3337
sales@macmedical.com
www.macmedical.com

MAC Medical, Inc.
Publication No. FFU-18 Rev D
July 2021
Information regarding this product is subject to change without prior notice.

