


Redesigned

Fixed and Adjustable Height Processing Sinks

Ergonomically
Friendly



Electronic adjustable height sink shown with optional lower shelves, pegboard backframe, and accessories

NEW FEATURES

- 12" Angled Backsplash
- Backsplash Mounted Faucets
- 30" Deep Standard Sink Top
- 20" Deep Standard Basin
- 2 1/2" Front of Sink to Basin
- Individual Plumbing for Adjustable Height Sinks
- Universal Workflow
- Larger Seismic Brackets
- Optional Pegboard Backframe and Accessories
- Optional Bolt-on Side Splash
- Optional H-Brace with Lower Shelf



Note: ETL compliance for adjustable height models only

[View more information on the back >>](#)

MADE
IN THE **USA**

mac
MEDICAL

STANDARD FEATURES

- Open Pipe Base
- Single, Dual, or Triple Basins
- Fully Welded Construction with Reinforced Tops
- 37" Standard Height for Stationary Sinks
- Adjustable Feet
- Marine Edge
- Flat Working Surface
- Etched Gallon Markings
- Basin Size: 24"W X 20"L X 10"D

OPTIONAL FEATURES

- **Ergonomically Friendly** Electronic Adjustable Height with Four Position Programming
 - Electronic Adjustable Height Range: 35"-42"
- Mixing Valve Faucet
- Pre-Rinse Faucet
- Deionized (DI) Water Faucet
- Air Gun Kit (with 5 tips)
- Water Gun Kit (with 8 tips)
- Eye Wash
- Larger Basin: 30"W X 20"L X 10"D
- Hand Lever Drains
- Copper Drain System (not compatible with deionized water)
- Back Frame

Back Frame Accessories (Sold Separately)

- Bin Rail
 - 3" Plastic Bin
 - 5" Plastic Bin
- Overhead LED Light
- Glove Dispenser
- Gallon Jug Holder
- Spray Bottle Holder
- Overhead Shelf
- Dosing System
- Monitor Mount

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at sales@macmedical.com.

Fixed Height Processing Sinks

Base Model #	Dimensions
PS307237-1000	30"D x 72"W x 37"H, Single Basin
PS309637-1000	30"D x 96"W x 37"H, Dual Basin
PS3012037-1000	30"D x 120"W x 37"H, Triple Basin

Note: Model # will change based on features added

Electronic Adjustable Height Processing Sinks

Base Model #	Dimensions
PS307235-1000	30"D x 72"W x 35"H, Single Basin
PS309635-1000	30"D x 96"W x 35"H, Dual Basin
PS3012035-1000	30"D x 120"W x 35"H, Triple Basin

Power Requirements (Electric Adjustable Height Sinks)

120 Volt, 60 Hz, single phase, 15 amp GFCI protected electrical outlet, installed per local building codes, and protective earth.

Note: Model # will change based on features added

INSPIRED
TO BE BETTER

Customer needs are our first priority.

MAC Medical, Inc.
Publication No. FFU-23 Rev C
April 2023

Information regarding this product is subject to change without prior notice.



Corporate Office
200 Admiral Weinel Blvd.
Columbia, Illinois 62236

Manufacturing Plant
820 South Mulberry Street
Millstadt, IL 62260

Phone: 618-476-3550 • Toll Free: 877-828-9975 • www.macmedical.com