


Fixed and Adjustable Height Processing Sinks

Ergonomically
Friendly



Electronic adjustable height sink shown
with optional lower shelves, pegboard
backframe, and accessories

STANDARD FEATURES

- 12" Angled Backsplash
- 37" Standard Height for Fixed Height Sinks
- Adjustable Feet
- Basin Size: 24"W X 20"L X 10"D
- Etched Gallon Markings
- Flat Working Surface
- Fully Welded Construction with Reinforced Tops
- Marine Edge
- Open Pipe Base
- Seismic Brackets
- Single, Dual, or Triple Basins



Note: ETL compliance
for adjustable height
models only

[View more information on the back >>](#)

MADE
IN THE **USA**

mac
MEDICAL

OPTIONAL FEATURES

(Must be specified with sink purchase)

- **Ergonomically friendly** electronic adjustable height hydraulic system with four position programming (electronic adjustable height range: 35" - 42")
- Bolt-on Side Splash
- Copper Drain System (Not Compatible with Deionized Water)
- Hand Lever Drains
- Larger Basin: 30"W X 20"L X 10"D

Fixed Height Processing Sinks

Base Model #	Dimensions
PS307237-1000	30"D x 72"W x 37"H, Single Basin
PS309637-1000	30"D x 96"W x 37"H, Dual Basin
PS3012037-1000	30"D x 120"W x 37"H, Triple Basin

Note: Model # will change based on features added

Electronic Adjustable Height Processing Sinks

Base Model #	Dimensions
PS307235-1000	30"D x 72"W x 35"H, Single Basin
PS309635-1000	30"D x 96"W x 35"H, Dual Basin
PS3012035-1000	30"D x 120"W x 35"H, Triple Basin

Power Requirements (Electric Adjustable Height Sinks)

120 Volt, 60 Hz, single phase, 15 amp GFCI protected electrical outlet, installed per local building codes, and protective earth.

Note: Model # will change based on features added

OPTIONAL SINK ACCESSORIES

- Backsplash Mounted Accessories
 - Air Gun Kit (with 5 Tips)
 - Deionized (DI) Water Faucet
 - Eye Wash
 - Gooseneck Mixing Valve Faucet
 - Pre-Rinse Faucet
 - Water Gun Kit (with 8 Tips)
- Basin Covers
- Leg Mounted Power Strip
- Lower Shelf/Shelves
- Wrist Rest

PEGBOARD BACKFRAME AND ACCESSORIES

- Pegboard Backframe
 - Bin Rail
 - Plastic Storage Bins – 2 Sizes Available
 - Brush Holder
 - Digital Basin Temperature Indicator
 - Dosing Soap Pump
 - Flushing System
 - Gallon Jug Holder
 - Glove Dispenser
 - Keyboard Arm
 - LED Magnifying Inspection Lamp
 - LED Soak Timer
 - Monitor Mount Arm
 - Overhead LED Light
 - Overhead Light Deflector
 - Pegboard Hooks
 - Pegboard Shelf
 - Spray Bottle Holder
 - Waterproof LED Task Light with Magnifying Lens
 - Work Flow Signs

To place an order, contact our customer service department at **618-476-3550** or **877-828-9975**, or by email at **sales@macmedical.com**.

INSPIRED
TO BE BETTER