

# TS-Series Touch Screen Blanket & Fluid Warming Cabinets




## Model # DWC183074T-TS

(T = Chambers in Thirds)


### DIMENSIONS

SHOWN WITH  
STANDARD 4" BASE




DOOR REMOVED  
FOR CLARITY

FRONT VIEW



SIDE VIEW



### SPECIFICATIONS

#### Overall Size

20.5"D x 30"W x 74.5"H  
(520.7mm x 762mm x 1892mm)  
12.27 Cubic Feet Storage  
3 Adjustable Shelves

#### Upper Chamber I.D.

17"D x 26"W x 15.25"H  
(431.8mm x 660.4mm x 387.3mm)  
3.45 Cubic Feet Storage  
Temperature Range 90°-160° F (Tolerance +/- 2° F)

#### Lower Chamber I.D.

17"D x 26"W x 34.5"H  
(431.8mm x 660.4mm x 876.3mm)  
8.82 Cubic Feet Storage  
Temperature Range 90°-160° F (Tolerance +/- 2° F)

#### Electrical

120V, 8.3 Amp, 50/60 Hz, .82 kWh (Avg.),  
2047 BTU/hr (Avg.)

#### Door Types

Stainless Steel or Glass