

Accessories

Magnifying Lamp



#H0509-09

LED Magnifying Inspection Lamp
with 12 and 5 Diopter lens



SPECIFICATIONS

Lens Size

12 Diopter Lens: 3.75 inch round
5 Diopter Lens: 4.75 inch round

Lens Material - Glass lens, rubber ring

Lens Diopter - 12 Diopter | 5 Diopter

Magnification

12 Diopter Lens: 4x magnification
5 Diopter Lens: 2.25x Magnification

Working Distance

12 Diopter Lens: Approx. 1 inch
5 Diopter Lens: Approx. 8 Inches

Arm Length (Reach) - 36 inch reach (max)

Arm and Lamp Adjustment - via tension knob

Color - Black

of LEDs - 64

LED Color Temperature - 6500K

Dimmable - YES

Expected LED life - Minimum 20,000 Hours

Lumen - Max 950 lm

Clamp Opening - 2.5" max

Power (Voltage)

Input: 100-240V~ 0.4A Max
Output: 26.0V --- 0.46A

LED Wattage - 1-10W

Power Cord Length - 5 Foot