


Accessories

Wrist Rests



#PS17E-24WR & PS17E-30WR

Ergonomic Wrist Rests for Processing Sinks



SPECIFICATIONS

PS17E-24WR

Fits 24" wide processing sink basin

Removable for easy cleaning

Size - 23.75 in W x 2.75 in H x 4.8 in D

Weight - 3.9 lbs

Material - Stainless steel

PS17E-30WR

Fits 30" wide processing sink basin

Removable for easy cleaning

Size - 29.75 in W x 2.75 in H x 4.8 in D

Weight - 4.9 lbs

Material - Stainless steel