

# Accessories

## Dosing Soap Pump



### Kit #PS17C-05

DoseTech DMS Automated Detergent Dosing System



Dosing Pump



Custom Bracket

## SPECIFICATIONS

**Brand** - KNIGHT

**Overall Dimensions** - 6" H x 7" W x 6-1/4" D  
(15.2 cm x 17.8 cm x 15.9 cm)

**Case Rating** - IP-64

**Dosing Pump Flow Rate** - 250 ml/min.

**Chemical Compatibility, Dosing Pump** -  
Industry standard enzymatic and other  
detergents

**Weight** - 5 lbs. (2.3 kg)

**Temperature Sensor** - Max temperature =  
257°F, 125°C

**Approvals** - CAN/CSA C22.2 No. 61010-1 ~  
UL61010-1

**Power Supply/Voltage** -

In: 110-240 VAC 1A 50-60 Hz

Out: 24 VDC 1.67A

