


# TS-Series Touch Screen Blanket & Fluid Warming Cabinets



Model # SWC182474-TS


## DIMENSIONS

SHOWN WITH STANDARD 4" BASE




DOOR REMOVED FOR CLARITY

FRONT VIEW



SIDE VIEW



## SPECIFICATIONS

### Overall Size

20.5"D x 24"W x 75"H  
(520.7mm x 609.6mm x 1905mm)  
2 Adjustable Shelves, 1 Fixed Shelf

### Chamber I.D.

17"D x 20"W x 61"H  
(431.8mm x 508mm x 1549mm)  
12 Cubic Feet Storage  
Temperature Range 90°-160° F (Tolerance +/- 3° F)

### Electrical

120V, 6.5 Amp, 50/60 Hz, .47 kWh (Avg.),  
1604 BTU/hr (Avg.)

### Door Types

Stainless Steel or Glass