


D-Series Data Logging Blanket & Fluid Warming Cabinets




Model # SWC182474

DIMENSIONS

SHOWN WITH STANDARD 4" BASE



DOOR REMOVED FOR CLARITY



SPECIFICATIONS

Overall Size

20.5"D x 24"W x 75"H
(520.7mm x 609.6mm x 1905mm)
2 Adjustable Shelves, 1 Fixed Shelf

Chamber I.D.

17"D x 20"W x 61"H
(431.8mm x 508mm x 1549mm)
12 Cubic Feet Storage
Temperature Range 90°-160° F (Tolerance +/- 3° F)

Electrical

120V, 6.5 Amp, 50/60 Hz, .47 kWh (Avg.),
1604 BTU/hr (Avg.)

Door Types

Stainless Steel or Glass