


# TS-Series Touch Screen Blanket & Fluid Warming Cabinets



Model # SWC243074-TS


## DIMENSIONS

SHOWN WITH STANDARD 4" BASE




DOOR REMOVED FOR CLARITY

FRONT VIEW



SIDE VIEW



## SPECIFICATIONS

### Overall Size

26.5"D x 30"W x 74.5"H  
(673.1mm x 762mm x 1892mm)  
2 Adjustable Shelves, 1 Fixed Shelf

### Chamber I.D.

23"D x 26"W x 61"H  
(584.2mm x 660.4mm x 1549mm)  
21.1 Cubic Feet Storage  
Temperature Range 90°-160° F (Tolerance +/- 3° F)

### Electrical

120V, 6.5 Amp, 50/60 Hz, .47 kWh (Avg.),  
1604 BTU/hr (Avg.)

### Door Types

Stainless Steel or Glass